



“ I’ve used CPAP, had surgeries and tried an oral device and none of them worked well for me. **My sleep has definitely improved with Provent.** – Bill ”



“ Doctors say it has given them a new weapon in the battle against sleep apnea, and **many patients who struggled with CPAP call (Provent) a godsend.** – Anahad O’Connor, The New York Times ”

“ For over a year, my wife tried everything and nothing has worked. In the last three days, **Provent has solved her problem—it’s just unbelievable.** – Carlos ”

“ When I first had Provent prescribed, my doctor warned me that it takes a little getting used to. Breathing out felt awkward at first. It took two nights, but by the third night I didn’t have the awkward feeling. **My doctor was right—I tried it long enough to get used to it and it worked great!** I’m very happy with how Provent has improved my sleep. – Marsha ”

EPAP: How Provent Works

When you inhale Provent’s MicroValve opens, allowing for nearly unobstructed airflow.



When you exhale, the MicroValve closes, restricting airflow to create EPAP (Expiratory Positive Airway Pressure) and maintain pressure in the airway through the start of the next inhalation.

For more information about Provent Therapy, ask your doctor.

Visit ProventTherapy.com
Call 1-888-SLP-WELL (1-888-757-9355)

Provent Sleep Therapy, LLC
373 South Willow Street, Suite 306
Manchester, New Hampshire 03103, USA


© 2016 Provent Sleep Therapy, LLC. Provent and the Provent logo are registered trademarks or trademarks of Provent Sleep Therapy, LLC in the U.S.A. and other countries. The Provent device is covered by U.S. and Foreign Patents: pat. pending and pat. <http://www.proventtherapy.com/patents>

PST008 Rev B

PROVENT[®]
sleep apnea therapy

PORTABLE
POWERFUL





Obstructive Sleep Apnea (OSA) is a serious medical condition characterized by pauses in breathing or shallow breaths during sleep.

Untreated OSA can have a negative impact on quality of life including decreased energy levels, lack of motivation, depression, and daytime sleepiness. It may even lead to additional health problems like high blood pressure, heart attack, stroke, and diabetes.

Simple

Provent Sleep Apnea Therapy is a simple, effective treatment option for OSA. Provent's proprietary MicroValves use the power of your own breathing to create Expiratory Positive Airway Pressure (EPAP) which helps keep your airway open, and enable natural airflow, while you sleep.

Effective & FDA Cleared

Available by prescription only, Provent Therapy is the only product of its kind FDA cleared for the treatment of OSA. In clinical studies Provent was shown to effectively treat mild, moderate and severe Obstructive Sleep Apnea.



Ideal for Travel

Small, discreet and disposable Provent patches are convenient for use nightly at home and while traveling. If you are non-compliant with CPAP therapy, or if you just want a more convenient OSA treatment option for travel, Provent just might be the answer for you.

American Made

Provent is proudly made in the USA. It is available with reimbursement through selected insurance companies, via private pay from local and national distributors, and directly to veterans through the VA Medical System. Visit proventtherapy.com/locator to find prescribers and suppliers near you.

Latex Free

Comfortable and breathable Provent patches are hypoallergenic, easy to apply, and they contain no natural rubber latex.



CPAP Mask & Machine Free

Provent Sleep Apnea Therapy has been on the market since 2012 and, with more than 5,000,000 nights prescribed, Provent is the proven CPAP alternative that requires no electricity, no mask and no machine!

The Provent Starter Kit

Like any new therapy Provent Sleep Apnea Therapy may take some getting used to. A 3-Phase Provent Starter Kit is available to help you get used to sleeping with the device, by gradually increasing EPAP resistance, over the first several nights of therapy.

Provent Therapy is a prescription product that should only be used after consultation with a physician, so talk to your doctor today and see if it could be the right solution for you.

As with any OSA treatment, there are risks associated with device use and the device is not effective in all patients. For additional product and safety information, please refer to the instructions for use included with the device.